



HEADPHONE SPLITTER

- If your child has a **wireless communication device** (e.g. Mini Mic, Roger Touchscreen and Receiver), you can plug this directly into the device you are accessing Online Services on so the audio signal can be streamed to your child's hearing device/s
- You can do this by directly plugging the wireless communication device into the audio jack/stereo jack (the little hole sometimes marked with a picture of headphones) using the **3.5mm audio cable (aux cable)**
- However, this means the sound would only go to your child's hearing devices. In order for you to hear the sound too, you will need to use something called a '**headphone splitter**' which splits the signal so you can plug your own headphones in too and hear the sound: you can think of it like a double adapter you might have used for power points
- This headphone splitter can be purchased easily and cheaply – if you search 'headphone splitter' online, you will see some options
- Please see the below diagram to help you visualise what the splitter does
- Please contact your audiologist if you have any questions about how to set this up for your child

